


## LAUNCH AND IMPLEMENTATION TIMING 2011-2012




### PHASE 1 : JULY 2011

New proximal-end connector compatible with current closed system containers.

*Incompatible with:*

- IV ports



### PHASE 2 : JANUARY 2012

New distal port design on ULTRAPAK® closed systems only compatible with the new SpikeRight® PLUS proximal-end connector.

*Incompatible with:*

- IV ports and spikes
- Former SpikeRight® design

Product	Product Code	Case UPC (GTN)	QTY/C
COMPAT DUALFLO® pump administration set with SpikeRight® PLUS connector and 1000 mL water bag	12154512	100-43900-89663-9	30
COMPAT® pump administration set with SpikeRight® PLUS connector	12154513	100-43900-89755-1	30
Enteral Spike Adapter with SpikeRight® PLUS Connector	12163246	100-43900-10662-2	36

**CAUTION:** Federal (USA) law and similar laws throughout the world restrict this device to sale by or on the order of a licensed healthcare practitioner.